



Recommended Maximum Pacing Duration for Defibrillator Pads (In Hours)

20ms Pulse Width Example: Physio Control, Life-Pack, 10P, 12, 15 and 20

		Pulse Amplitude (mA)																			
		10	20	30	40	50	60	70	80	90	100	110	120	130	140	150	160	170	180	190	200
Pacing Rate (ppm)	40	24	24	24	24	24	24	24	22	20	18	18	16	14	12	11	10	8	8	8	8
	50	24	24	24	24	24	22	23	21	20	18	17	16	14	12	10	8	8	7	6	6
	60	24	24	24	24	22	20	20	20	18	17	16	15	13	11	8	6	4	4	5	4
	70	24	24	24	22	20	18	17	17	16	16	15	13	11	10	6	5	3	3	4	3
	80	24	24	24	22	18	16	15	14	14	13	12	11	10	9	6	4	3	3	3	3
	90	24	24	23	18	16	14	14	13	13	12	10	10	9	8	5	3	3	2	2	2
	100	24	24	20	18	16	14	12	12	11	10	9	8	8	7	5	3	2	2	2	2
	110	24	24	16	16	12	11	11	11	9	9	8	8	7	6	4	3	2	2	2	1
	120	24	24	16	14	12	10	9	8	8	7	6	5	5	4	4	2	2	2	1	1
	130	24	24	14	14	10	9	8	7	6	5	4	4	4	4	3	2	2	1	1	1
	140	24	24	12	10	10	8	6	5	4	4	4	3	3	3	2	2	2	1	1	1
	150	24	24	10	10	8	6	5	4	4	4	3	3	2	2	2	2	1	1	1	1
	160	24	24	8	8	8	5	5	4	3	3	3	2	2	2	2	1	1	1	1	1
170	24	24	8	8	6	5	4	4	3	3	3	2	2	2	2	1	1	1	1	1	



Recommended Maximum Pacing Duration for Defibrillator Pads (In Hours)

40ms Pulse Width Example: Zoll

		Pulse Amplitude (mA)																				
		10	20	30	40	50	60	70	80	90	100	110	120	130	140	150	160	170	180	190	200	
Pacing Rate (ppm)	30	24	24	24	22	20	16	14	12	10	9	8	8	7	6							
	40	24	24	24	16	15	14	12	10	9	8	7	6	5	5							
	50	24	24	22	14	14	12	10	9	8	7	6	6	5	4							
	60	24	24	20	12	12	10	10	9	8	6	6	5	5	4							
	70	24	24	16	10	10	9	8	8	6	5	4	4	3	3							
	80	24	24	14	8	8	7	6	6	4	4	3	3	2	2							
	90	24	24	12	8	7	6	5	5	4	3	3	2	2	2							
	100	24	24	10	6	6	5	5	4	4	3	3	2	2	2							
	110	24	24	8	6	5	4	4	3	3	3	2	2	2	2							
	120	24	24	8	7	4	4	3	3	3	2	2	2	1	1							
	130	24	24	7	6	4	4	3	3	2	2	2	2	1	1							
	140	24	24	6	5	3	3	3	2	2	2	1	1	1	1							
	150	24	24	6	4	3	3	2	2	2	2	1	1	1	1							
	160	24	24	5	4	3	2	2	2	2	2	1	1	1	1							
	170	24	24	5	4	2	2	2	2	1	1	1	1	1	1							
	180	24	24	4	3	2	2	2	2	1	1	1	1	1	1							



Recommended Maximum Pacing Duration for Defibrillator Pads (In Hours)

20ms Pulse Width Example: Philips Codemaster

Pulse Amplitude (mA)																					
	10	20	30	40	50	60	70	80	90	100	110	120	130	140	150	160	170	180	190	200	
Pacing Rate (ppm)	40	24	24	24	24	24	24	24	22	20	18	18	16	14	12	11	10	8	8		
	50	24	24	24	24	24	22	23	21	20	18	17	16	14	12	10	8	8	7		
	60	24	24	24	24	22	20	20	20	18	17	16	15	13	11	8	6	4	4		
	70	24	24	24	22	20	18	17	17	16	16	15	13	11	10	6	5	3	3		
	80	24	24	24	22	18	16	15	14	14	13	12	11	10	9	6	4	3	3		
	90	24	24	23	18	16	14	14	13	13	12	10	10	9	8	5	3	3	2		
	100	24	24	20	18	16	14	12	12	11	10	9	8	8	7	5	3	2	2		
	110	24	24	16	16	12	11	11	11	9	9	8	8	7	6	4	3	2	2		
	120	24	24	16	14	12	10	9	8	8	7	6	5	5	4	4	2	2	2		
	130	24	24	14	14	10	9	8	7	6	5	4	4	4	4	3	2	2	1		
	140	24	24	12	10	10	8	6	5	4	4	4	3	3	3	2	2	2	1		
	150	24	24	10	10	8	6	5	4	4	4	3	3	2	2	2	2	1	1		
160	24	24	8	8	8	5	5	4	3	3	3	2	2	2	2	1	1	1			
170	24	24	8	8	6	5	4	4	3	3	3	2	2	2	2	1	1	1			
180	24	24	8	8	6	5	4	3	2	2	2	1	1	1	1	1	1	1			



Recommended Maximum Pacing Duration for Defibrillator Pads (In Hours)

40ms Pulse Width Example: Philips HeartStart MRx Monitor/Defibrillator/Pacemaker

		Pulse Amplitude (mA)																				
		10	20	30	40	50	60	70	80	90	100	110	120	130	140	150	160	170	180	190	200	
Pacing Rate (ppm)	30	24	24	24	22	20	16	14	12	10	9	8	8	7	6	5	4					
	40	24	24	24	16	15	14	12	10	9	8	7	6	5	5	4	4					
	50	24	24	22	14	14	12	10	9	8	7	6	6	5	4	4	3					
	60	24	24	20	12	12	10	10	9	8	6	6	5	5	4	3	3					
	70	24	24	16	10	10	9	8	8	6	5	4	4	3	3	3	3					
	80	24	24	14	8	8	7	6	6	4	4	3	3	2	2	2	2					
	90	24	24	12	8	7	6	5	5	4	3	3	2	2	2	2	1					
	100	24	24	10	6	6	5	5	4	4	3	3	2	2	2	2	1					
	110	24	24	8	6	5	4	4	3	3	3	2	2	2	2	1	1					
	120	24	24	8	7	4	4	3	3	3	2	2	2	1	1	1	1					
	130	24	24	7	6	4	4	3	3	2	2	2	2	1	1	1	1					
	140	24	24	6	5	3	3	3	2	2	2	2	1	1	1	1	1	1				
	150	24	24	6	4	3	3	2	2	2	2	1	1	1	1	1	1	1				
160	24	24	5	4	3	2	2	2	2	2	1	1	1	1	1	1	1					
170	24	24	5	4	2	2	2	2	2	1	1	1	1	1	1	1	1					
180	24	24	4	3	2	2	2	2	2	1	1	1	1	1	1	1	1					