



Specialty Defibrillation Electrodes

10 Foot Lead Wire: Zoll




10' lead wire allows the extra distance the clinician needs in any room configuration or clinical situation.

Product: T100AC-Zoll-10

(10/box, \$240/box)

C200: Physio, Philips, Zoll



Our specialty, transparent electrode offers a larger back/posterior pad to bringing down impedance. Thus, decreasing the chance of patient burns.

**Products: C200AC-Physio
C200AC-Philips
C200AC-Zoll**

(10/box, \$380/box)

C300: Physio, Philips Zoll



This specialty, transparent electrode provides increased surface area that giving these electrodes the maximum amount of adhesion for diaphoretic patients.

**Products: C300AC-Physio
C300AC-Philips
C300AC-Zoll**

(10/box, \$240/box)