



HEARTSYNC®

Defibrillator & Electrode Compatibility

Brand of Defibrillator	Model of Defibrillator	Heart Sync Compatible Defibrillation Electrode Part Numbers
Philips Medical	» XLT » MRX » HeartStart XL » MRX ALS » HP Codemaster XL (Plug Style Connector)	» T100AC-Philips » C100LOAC-Philips » T100LOAC-Philips » Pediatric-Philips » C100AC-Philips
	» HeartStart FR2+ » HeartStart Forerunner	» T100AC-Philips » C100AC-Philips » T100LOAC-Philips » C100LOAC-Philips
	» HeartStart FR3 » HeartStart FRx	Heart Sync defibrillation electrodes are NOT compatible with these Philips devices
Physio-Control	» LP-9P » LP-15 » LP-10C » LP-20 » LP-12	» T100AC-Physio » C100LOAC-Physio » T100LOAC-Physio » Pediatric-Physio » C100AC-Physio
	» LP-500 » LP-500 DPS » LP-1000	» T100AC-Physio » C100AC-Physio » T100LOAC-Physio » C100LOAC-Physio
Zoll Medical	» M-Series » E-Series » R-Series » X-Series	» T100AC-Zoll » C100LOAC-Zoll » T100LOAC-Zoll » Pediatric-Zoll » C100AC-Zoll
	» Zoll AED PLUS » Zoll AED Pro	Heart Sync defibrillation electrodes are NOT compatible with Zoll AED's.
Cardiac Science	» PowerHeart AED G3 Plus » PowerHeart AED G3 Pro	» T100-CS

Heart Sync's pediatric defibrillation electrodes are NOT compatible with AED's. AED's cannot be manually adjusted to reduce energy delivery. Please refer to your AED's manufacturer operation manual.

Contact Heart Sync directly for any additional questions regarding electrode compatibility.