





Placement and User Recommendations

Confirm compatibility before use

Defibrillation Only




Pacing and Defibrillation Placement





Guidelines For Child Use:

Manual Defibrillator:
≥ 10 kg (22lbs)


AED or AED Mode:
≥ 25 kg (55lbs)
≥ 8 years old



Defibrillation Only




Pacing and Defibrillation Placement



Guidelines For Infant Use:

Manual Defibrillator:
≤ 10 kg (22lbs)

AED (with a red 'X' over the text)



ADULT/CHILD
≥ 10kg (22 lbs)

ADULT/CHILD
≥ 10kg (22 lbs)

INFANT (Pediatric)
< 10kg (22 lbs)

INFANT (Pediatric)
< 10kg (22 lbs)